MFP103B

GEAR PUMP

Gear pump for use with the Positive Displacement Pump Support Module (MFP103)





- Popular design pump for use with TecQuipment's Positive Displacement Pump Support Module (MFP103)
- Quick-release, self-sealing connections for simple and safe fitting
- Shows the characteristics of a gear pump



GEAR PUMP

DESCRIPTION

For use with the Positive Displacement Pump Module (MFP103) this pump is ideal for student experiments, demonstrations and projects.

The gear pump is a positive displacement pump. It has double gears that deliver a given volume of fluid (oil) for each full rotation of the pump shaft.

The gear rotation generally determines the flow direction, but you only test the pump in one direction, determined by the Universal Dynamometer.

Self-sealing connections reduce oil spillage and simplify installing the pump to the pump module.

STANDARD FEATURES

- Five-year warranty
- Made in accordance with the latest European Union directives

LEARNING OUTCOMES

- Performance and characteristics of a gear pump
- · Volumetric and overall efficiencies
- Use of an oval gear flowmeter

WHEN TWO OR MORE OPTIONAL PUMPS ARE ORDERED:

 Comparison of positive displacement pumps (economy, flow rate and output pressure pulses)

ESSENTIAL BASE UNIT

 Positive Displacement Pump Module (MFP103) (with Universal Dynamometer MFP100)

OPERATING CONDITIONS

OPERATING ENVIRONMENT:

Laboratory

STORAGE TEMPERATURE RANGE:

-25°C to +55°C (when packed for transport)

OPERATING TEMPERATURE RANGE:

+5°C to +40°C

OPERATING RELATIVE HUMIDITY RANGE:

80% at temperatures < 31°C decreasing linearly to 50% at 40°C

SPECIFICATION

NETT DIMENSIONS:

Approximately 250 mm x 250 mm x 250 mm and 8 kg

PACKED DIMENSIONS:

Approximately 0.05 m³



TECQUIPMENT LTD, BONSALL STREET, LONG EATON, NOTTINGHAM NGIO 2AN, UK TECQUIPMENT.COM +44 115 972 2611 SALES@TECQUIPMENT.COM

DB 0517